


CAESA-LAB INC.  
www.caesalabinc.com

# BACTRON CHAMBERS



CAESA-LAB INC.  
info@caesalabinc.com  
www.caesalabinc.com  
Tel: (450) 441-5408

## BACTRON Anaerobic/Environmental Chambers

A Cost Saving Solution for Improved Anaerobe Culture & Identification



Patented BACTRON anaerobic/environmental systems allow efficient and dexterous glove-free handling and inspection of samples. Comfortable arm cuffs, which are available in three sizes, fit snugly around the users forearm. The convenient foot pedals are used to evacuate oxygen and add anaerobic gas to the sleeves while the airlock auto-cycle facilitates oxygen-

free transfer of samples and supplies in the chamber. Indicator lights will signal that the air-lock is anaerobic and ready to be opened into the chamber. Simply open arm port doors and place on the convenient shelf or racks for out-of-the-way storage. Samples and supplies can then be brought into the chamber with ease from the sliding airlock shelf.

### Cost Comparison: Anaerobic Jars vs. Bactron III

The bottom line in operating costs:



Anaerobic jars cost \$0.27 per day, per plate.

The Bactron III costs only \$0.05 per day per plate to operate



300 Plates/Day	Anaerobic Jars		Bactron III Anaerobic Chamber	
	Annual Cost	Five Year	Annual Cost	Five Year
Consumables:				
30 Jars per day @ \$2.46 (gas generating kit, strips catalyst sachets)	\$26,937.00	\$134,685.00		
Gas, test strips, sleeves, catalyst			\$2,007.50	\$10,037.50
Equipment:				
Jars				
30 - 10 Plate Jars @ \$312.43 ea.	\$9,372.90	\$9,372.90		
2 Incubators	\$6,500.00	\$6,500.00		
Bactron III Anaerobic Chamber			\$16,000.00	\$16,000.00
Total:		\$150,557.90		\$26,037.50
Five Year Equipment and Operating Costs				
Per Day Cost Per Plate		\$0.27		\$0.05

# BACTRON Anaerobic/Environmental Chambers

## An Overview of the Unique BACTRON Features



Vacuum pump, hoses and fittings are included.

The air-lock feature offers ample pass-through capacity with a sliding shelf (shelf comes with II, III, IV only)

Versatile, glove-free cuff system allows dexterous handling for users of any size.

Optional powder-coated, stand with locking casters.

Durable, powder-coated, stainless steel body with rigid plexiglas front for clear viewing of samples.

Internal Incubator Ranging from 100 - 900 plates depending on the model.

Large, ambient temperature workspace for user comfort and convenience.

### Features/Benefits:

- Gloveless design for dexterous, easy handling of samples
- Separate incubator and work chamber for increased efficiency and user comfort
- Condensate controller for virtually condensation-free work area and fog-free viewing
- Catalyst cartridge increases anaerobic results and efficiency
- System ready for out-of-the-box use, just add gas
- Airlock can be operated manually or via unique auto-cycle
- Microprocessor temperature controller with accurate and user-friendly keypad calibration
- Independent over-temperature safety controller
- Easy-to-use sliding airlock shelf preserves work space
- Chamber pressure relief system featuring BACTRON manometer for comfort and maintenance of anaerobic environment

### Included with each system:

- Vacuum pump with fittings and hoses
- One year warranty\*
- Two Catalyst Cartridges
- Petri plate racks
- Sleeve assembly with medium cuff (large and small cuffs are also available)

### Education:

Each BACTRON purchase within the USA, includes set up and on-site training by the local SHEL LAB Sales Representative. One staff member is eligible to attend a seminar on Anaerobe identification procedures and techniques.\*\*



# BACTRON Anaerobic/Environmental Chambers

## An Overview of the Unique BACTRON Features

### BACTRON: Value

The BACTRON Anaerobic Chamber systems are just that: a system. All components for trouble-free use are included with your initial purchase, including a vacuum pump, sleeves, cuffs and all the necessary tubing.


### BACTRON: Versatility

offers multiple high-value features and essential benefits. Incubator Petri plate capacities and chamber dimensions accommodate varied laboratory needs and budgets. Ranging from 100 to 900 Petri plate incubator capacities, the BACTRON Chamber line can address a myriad of anaerobic, hypoxic or other special gas environment applications.

### BACTRON: Innovation

The dedicated incubators mean that all Bactron Chambers have an ambient-temperature workspace. This allows users to stay inside the chamber for extended periods of time without becoming uncomfortable. Our innovative humidity control system provides a condensation free environment that allows for clear, easy inspection of cultures.

The Bactron Chamber is constructed of highly resilient materials. The sturdy Plexiglas hood and durable stainless steel body ensure years of long, trouble-free use. The BACTRON line of anaerobic chambers feature a unique pressure-management system inside the chamber for comfortable use and strict anaerobic conditions. A manometer is located inside the chamber to provide a visual reference and to relieve pressure during use for comfort and convenience.



### Quick Tip



Many BACTRON users have learned to save time, and conserve gas by bringing small items into the chamber using the gloveless sleeve system.

### Airlock Tray



The Interior Airlock of the BAC II, III, and IV is equipped with a sliding tray to easily bring items into the anaerobic workspace.

### Applications:

- Anaerobic Microbiology Assays
- Clinical Microbiology
- Anaerobic Tissue Culture Research
- Anaerobic Bacteriology
- Hypoxic Tissue Culture Research
- Microaerophilic/Hypoxic
- Quality Control
- Packaging
- R&D
- Industrial Research

BACTRON Model <i>*-2 denotes 230V</i>	BACTRON I BACTRON I-2	BACTRON II BACTRON II-2	BACTRON III BACTRON III-2	BACTRON IV BACTRON IV-2	BACTRON X BACTRON X-2
Exterior Dimensions WxDxH	48.8" x 32.3" x 27.5" 124 x 82 x 69.9 cm	61.3" x 32.3" x 27.5" 155.7 x 82 x 69.9 cm	48.8" x 32.3" x 33.5" 124 x 82 x 85.1 cm	87.8" x 32.3" x 27.5" 223 x 82 x 69.9 cm	61.3" x 32.3" x 27.5" 155.7 x 82 x 69.9 cm
Volume	43346.6 / 710.3L	54449.7 / 892.3L	52804 / 865.3L	77988.4 / 1278L	54444.9 / 892.3L
Interior Dimensions	33" x 28.5" x 25"* 83.8 x 72.4 x 63.5 cm	42.3" x 28.5" x 25"* 107.4 x 72.4 x 63.5 cm	33" x 28.5" x 25"* 83.8 x 72.4 x 63.5 cm	42.3" x 28.5" x 25"* 107.4 x 72.4 x 63.5 cm	42.3" x 28.5" x 25"* 107.4 x 72.4 x 63.5 cm
Volume	21136.5 / 346.4L	27093.2 / 444L	21136.5 / 346.4L	27093.2 / 444L	27093.2 / 444L
Air Lock (Access)	9" x 11" x 9" 22.9 x 27.9 x 22.9 cm	12" x 14" x 12" 30.5 x 35.6 x 30.5 cm	9" x 11" x 9" 22.9 x 27.9 x 22.9 cm	12" x 14" x 12" 30.5 x 35.6 x 30.5 cm	12" x 14" x 12" 30.5 x 35.6 x 30.5 cm
Volume	1053 / 17.3L	2304 / 37.8L	1053 / 17.3L	2304 / 37.8L	2304 / 37.8L
Incubator WxHxD	26.8" x 9.25" x 8.3" 68.1 x 23.5 x 21.1 cm	36" x 9.25" x 8.3" 91.4 x 23.5 x 21.1 cm	27.5" x 13.6" x 8.5" 69.9 x 34.5 x 21.6 cm	23" dia x - x 18.5" 58.4 dia x - x 47 cm	N/A
Volume	2057.6 / 33.7L	2763.9 / 45.3L	3179 / 52.1L	7686.3 / 126L	
<b>Electrical</b>					
Voltz/Hertz	120-230/50-60	120-230/50-60	120-230/50-60	120-230/50-60	120-230/50-60
Watts	1100	1100	1100	1350	1100
Phase/amps	Single/9.17	Single/9.17	Single/9.17	Single/11.25	Single/9.17
Interior Outlets	2** 1amp. Max	2** 1amp. Max	2** 1amp. Max	2** 1amp. Max	2** 1amp. Max
Exterior Outlet, vacuum	5 amp. Max	5 amp. Max	5 amp. Max	5 amp. Max	5 amp. Max
Circuit breakers main/vacuum	2	2	2	2	2
cUL, CE Approved	Yes	Yes	CE Only	Yes	N/A
<b>Incubator Performance</b>					
Temp. range	Amb +5° to 70°C	Amb +5° to 70°C	Amb +5° to 70°C	Amb +5° to 70°C	N/A
Temp. uniformity	(+/-) .25° C	(+/-) .25° C	(+/-) .25° C	(+/-) .25° C	N/A
<b>Weight</b>					
Shipping Weight lb (kg)	480 (218)	610 (277)	480 (218)	760 (345)	610 (277)
With optional stand lb (kg)	710 (323)	760 (346)	710 (323)	1310 (596)	760 (346)
<b>Miscellaneous</b>					
Petri plate racks included	6	8	7	10	N/A
Optional stand dimensions	49" x 30" x 30" 124.5 x 76.2 x 78.7 cm	61.5" x 30" x 30" 156 x 76.2 x 78.7 cm	49" x 30" x 30" 124.5 x 76.2 x 78.7 cm	88.5" x 30" x 31" 225 x 76.2 x 78.7 cm	61.5" x 30" x 30" 156 x 76 x 78.7 cm
Plate capacities	100	200	300	600 (900 with upgrade)	N/A

Catalyst Oven



Optional Accessories

# BACTRON Anaerobic/Environmental Chambers

## The Bactron Product Offering

### BACTRON I



The BACTRON I & II are very similar units differing in size but not in performance. The BACTRON I is the smallest of the fleet with a 100 plate capacity. Our most popular seller, the BACTRON II, provides larger workspace and plate capacity.

### BACTRON II



The BACTRON III uses a rear wall incubator with sliding doors to increase plate capacity to 300. This compact design provides the most incubation room with the smallest foot print

### BACTRON III



The BACTRON IV offers 600 plate capacity by using a space saving rotating shelf design. Two rear mounted shelves help keep your workspace organized.

### BACTRON IV



### BACTRON IV 900



The BACTRON IV 900 has the same specifications as the standard BACTRON IV with the addition of a rear wall incubator. With capacity of 900 plates this unit will suffice for even the busiest lab.