

SPECIMEN COLLECTION

Specimen Collection Kit

For specimen collection and transport, this system provides disposable products for the convenient collection, delivery, and testing of urine specimens. A plastic cup is provided for specimen collection and the pressure sensitive patient identification labels have ample writing space. The lightweight transport rack is small and easy to carry yet provides space for up to 20 tubes with labels and caps. The tubes are made of sturdy polystyrene and are safe for centrifuging. The flared top permits the use of a midget urinometer ensuring a tight cap fit. All elements of the system are disposable. A case contains 500 specimen collection cups, 500 15 mL collection tubes, 500 white caps, 500 pressure sensitive labels, and 5 transport racks.



Specimen Collection Kit

24 Hour Urine Container BMP-UR-554

This polyethylene container features easy-to-read graduations, a convenient carrying handle and an amber color for UV-protection. The plastic cap has an extra-thick F-217 liner to prevent leakage. The containers are easily stored and take minimal space in the refrigerator. Inert to normal preservatives. Container volume is 3000 mL.

24 Hour Urine Container with Dripless Pour Spout BMP-UR-552

This 3000 mL polyethylene container assures the precise measuring and dispensing of urine samples. The precision-made, fine pour spout allows for the sanitary transfer of small sample quantities into secondary containers. An integral dripless cap prevents cross-contamination and is attached to prevent cap loss. The container features easy-to-read graduations, a convenient carrying handle, and an amber color for UV protection.



Sputum Collection System



BMP-UR-552

BMP-UR-554

Sputum Collection System

Our sterile system is preassembled and packaged in an easy opening polyethylene bag. All components are polypropylene and are fully autoclavable. Our system consists of a wide funnel, lid, screw cap, graduated 50 mL centrifuge tube, identification label, and a base unit.



Sputum Collection Packaging

Item Number	Description	Case Pack
URINALYSIS SYSTEM		
BMP-UR-570	Specimen Collection Kit	500
TUBES AND CAPS IN BULK		
BMP-UR-572	Specimen Urinalysis Tubes only	500
BMP-UR-576	Caps only	1000
24 HOUR URINE CONTAINERS		
BMP-UR-554	24 Hour Urine Container	40
BMP-UR-552	24 Hour Urine Container with Dripless Spout	40
SPUTUM COLLECTION KIT		
BMP-SP-102	Sputum Collection Kit	36

CAESA-LAB INC.

Tel: (450) 441-5408
www.caesalabinc.com

SPECIMEN COLLECTION



Commode Specimen Collection System

This system makes stool collection easier and cleaner. It has easy-to-read graduations and a closed container transportation system. Its sturdy frame enables the container to fit securely under any toilet seat. The container has a removable, support frame to allow for easier transport to the lab. The directions are printed on the lid of the container in both English and Spanish. This system may also be used for female urine collection.



Sterile Midstream Kit

Sterile Midstream Kits

The 5 oz. widemouth polypropylene sterile sample container has a leakproof Plastisol lined metal cap and is provided with a patient identification label and towelettes. Fully autoclavable.



Sterile Wide Mouth Specimen Containers

Sterile Wide Mouth Specimen Containers

This convenient 5 oz. size allows for easy handling. A metal cap with Plastisol liner provides a leakproof seal.

Item Number	Approx. Capacity		Case Pack
COMMODOE SPECIMEN COLLECTION SYSTEM			
BMP-CN-894	22 oz.	700 cc	60
STERILE MIDSTREAM KITS			
BMP-UR-550	5 oz.	147.87 mL	100
STERILE WIDE MOUTH SPECIMEN CONTAINERS			
BMP-CN-890	5 oz.	147.87 mL	200

CAESA-LAB INC.

Tel: (450) 441-5408
www.caesalabinc.com