

NEW

CAESA-LAB

Mini Centrifuge



- + Eight position microtube rotor
- + Storage compartment for PCR rotor
- + Near silent operation
- + Starts and stops in just seconds
- + Conserves valuable bench space

Technical Data

Speed:	6,000rpm / 2,000 xg
Capacity:	8 x 1.5 / 2.0ml tubes 16 x 0.2ml PCR tubes 2 x 0.2ml PCR strips
Dimensions:	4.5 x 5.9 x 4.5 in 11.4 x 15 x 11.4 cm
Weight:	1.1 lbs / 5 kg
Electrical:	100 - 240V, 50-60Hz

Eight Tube Capacity

The new CAESA-LAB Mini centrifuge is a personal centrifuge that virtually fits in the palm of your hand. Ideal for quick spin downs of microtubes and PCR tubes, the centrifuge is extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpm. Open the lid, and the rotor quickly decelerates to a stop for removing samples.

An increased capacity eight position microtube rotor is included and comes installed. A second rotor for 0.2ml PCR tubes and strips is also included and stored with the rotor removal key in a convenient compartment on the bottom of the centrifuge. When the rotor is not in use, storage in this compartment prevents loss or misplacement of the rotor.



Storage Compartment



Includes Both Rotors

Ordering Information

CA-1008-B	Mini Centrifuge Blue lid
CA-1008-C	Mini centrifuge Clear lid
CA-1008-G	Mini centrifuge Neon Green lid
CA1-008-P	Mini centrifuge Purple lid
CA-1008-R	Mini centrifuge Red lid

*All include microtube rotor, PCR rotor, 0.5ml & 0.2ml adapters

CAESA-LAB

Tel: (450) 441-5408

www.caesalabinc.com