


## EMBEDDING TISSUE CASSETTES



CAESA-LAB disposable plastic cassettes can hold tissue specimens during the embedding process, as well as in a storage file. Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. They have a writing anterior area at a 30° angle.

Item No.	Specification	Quantity (pcs)	Measurement (cm)
F1009W	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009Y	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009G	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009P	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009LB	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009DB	Tissue Embedding Cassette	5000	48 X 32 X 44
F1010	Embedding Cassette	4000	48 X 32 X 32

### CAESA-LAB INC.

TELEPHONE : (450) 441-5408  
 FAX : (450) 441-9167  
 E-MAIL : info@caesalabinc.com  
 www.caesalabinc.com