

# SM-Plus Series

Dual Magnification Deluxe Stereo Microscope

# Microscope



## Features

- Ergonomic handle makes transporting microscope safe and secure
- Safe and durable dove-tail stand
- Super bright illumination with 30 LED bulbs for incident light
- Unique stage plate creates uniform transmitted illumination
- Adjustable top and bottom illumination
- Extra long working distance

## Specifications



Model	SM-Plus-13	SM-Plus-24
<b>Specification</b>		
Head	45° inclined binocular head	
Eyepiece	WF 10X/20	
Working Distance	100mm	
Interpupillary Distance	55-75mm	
Illumination	Reflected (30LED Bulbs) and Transmitted LED illumination	
Power	85V-265V	
Objectives	1X and 3X	2X and 4X
Total Magnification	10X and 30X	20X and 40X

Objective Magnification	Visual Field
1X	20.0 mm
2X	10.0 mm
3X	7.5 mm
4X	5.0 mm