

CAESA-LAB INC.


Tel: (450) 441-5408 / Fax: (450) 441-9167
www.caesalabinc.com e-mail: info@caesalabinc.com



Plastic Vials			Metric		Standard	
Part No. (Dram)	Fluid Ounce	Cubic cm (ml's)	Inside Diameter	Inside Height	Inside Diameter	Inside Height
2.5U	.28 oz	8.4 ml	17.0 mm	36.5 mm	10/16"	1 7/16"
3U	.40 oz	11.80 ml	17.0 mm	52.4 mm	10/16"	2 1/16"
3.5U	.44 oz	12.94 ml	20.6 mm	36.5 mm	13/16"	1 7/16"
5C	.63	18.48	20.6	52.4	13/16"	2 1/16"
7C	.88	25.88	24.9	52.4	1"	2 1/16"
9C	1.13	33.27	24.9	66.7	1"	2 5/8"
12C	1.50	44.36	28.6	66.7	1 1/8"	2 5/8"
15C	1.88	55.45	28.6	82.6	1 1/8"	3 1/4"
20C	2.50	73.93	38.1	68.3	1 1/2"	2 11/16"
25C	3.13	92.42	38.1	84.1	1 1/2"	3 5/16"
40C	5.00	147.86	47.6	84.1	1 7/8"	3 5/16"
50C	6.25	184.83	47.6	108.0	1 7/8"	4 1/4"

These transparent crystal polystyrene vials with tabbed white polyethylene snap-caps include 12 sizes from .28 ounces to 6.25 ounces.

These vials are ideal for use in the following industries: Pharmaceutical, Medical, Agricultural / Toxicology, Hobbyist, Home Organization.



Vials are shown as actual size. Vial dimensions are approximate Made in the U.S.A.