


Racks

RACKS FOR TEST TUBES

Our well-designed test tube racks are made of high density polypropylene, and are autoclavable and non-floating in a water bath. Sizes range from 13mm to 30mm and, therefore, accommodate all routine tubes. The 30mm size rack is especially suitable for 50mL centrifuge tubes. Alphanumeric markings are molded into the racks to provide easy test tube identification. They are available in four colors for color coding of test samples, or simply for aesthetics in the modern laboratory. The racks are shipped flat in order to reduce shelf space and then simply snapped together for use. They are stackable for convenient storage.

- Produced from polypropylene
- Autoclavable at 12 1°C for 20 minutes
- Non-floating in a water bath.
- Can be used in freezers without deforming or becoming brittle



For 13mm Tubes

90 holes, 6 x 15 format, rack dimensions: 246 x 105 x 64mm

| Item # | Description | Unit |
|--------|-----------------------------|------|
| 456403 | 13mm size, 90-place, white | Each |
| 456404 | 13mm size, 90-place, blue | Each |
| 456406 | 13mm size, 90-place, yellow | Each |
| 456410 | 13mm size, 90-place, red | Each |

For 16mm Tubes

60 holes, 5 x 12 format, rack dimensions: 246 x 105 x 72mm

For 25mm Tubes

| Item # | Description | Unit |
|--------|-----------------------------|------|
| 456503 | 16mm size, 60-place, white | Each |
| 456504 | 16mm size, 60-place, blue | Each |
| 456506 | 16mm size, 60-place, yellow | Each |
| 456510 | 16mm size, 60-place, red | Each |

For 20mm Tubes

40 holes, 4 x 10 format, rack dimensions: 246 x 105 x 72mm

| Item # | Description | Unit |
|--------|-----------------------------|------|
| 456603 | 20mm size, 40-place, white | Each |
| 456604 | 20mm size, 40-place, blue | Each |
| 456606 | 20mm size, 40-place, yellow | Each |
| 456610 | 20mm size, 40-place, red | Each |

40 holes, 4 x 10 format, rack dimensions: 295 x 125 x 85mm

| Item # | Description | Unit |
|--------|-----------------------------|------|
| 456703 | 25mm size, 40-place, white | Each |
| 456704 | 25mm size, 40-place, blue | Each |
| 456706 | 25mm size, 40-place, yellow | Each |
| 456710 | 25mm size, 40-place, red | Each |

For 30mm Tubes

24 holes, 3 x 8 format, rack dimensions: 300 x 112 x 85mm

| Item # | Description | Unit |
|--------|-----------------------------|------|
| 456803 | 30mm size, 24-place, white | Each |
| 456804 | 30mm size, 24-place, blue | Each |
| 456806 | 30mm size, 24-place, yellow | Each |
| 456810 | 30mm size, 24-place, red | Each |

RACKS FOR MICROCENTRIFUGE TUBES

Caesa's sturdy microcentrifuge tube racks are made of high density polypropylene and hold up to 100 microcentrifuge tubes. They are autoclavable and are non-floating in a water bath. Alphanumeric grid marks are molded into the racks to identify tube positions. The racks are supplied flat to reduce shelf space and transportation costs. Racks are stackable and have spaces for labels on each side (for bar coding, etc...)

- 100 holes, rack dimensions: 263 x 109 x 45mm

| Item # | Description | Unit |
|--------|---|------|
| 456303 | Microcentrifuge Tube Rack, 100-place, white | Each |

