

Bags

SPECIMEN TRANSPORT BAGS - ZIPPER CLOSURE

CAESA-LAB'S specimen transport bags are designed for the safe transportation of specimens and documentation. Each bag is printed with the biohazard symbol and offers checkmark boxes to indicate special handling instructions.

- 2 mil gauge polyethylene (PE)
- Attached document compartment
- Zipper closure allows repeated access to specimen
- Imprinted with biohazard symbol and score line
- Complies with OSHA and NCCLS guidelines
- Item #4922 features an absorbent pad that can hold 60mL of urine and 35mL of blood.



Item #	Description	Unit
4918	6 x 9" specimen transport bag, zipper closure	1000
4922	6 x 9" specimen transport bag, zipper closure, with absorbent pad	500
4950A	6 x 9" STAT specimen transport bag, zipper closure, red tint	1000
4926	8 x 10" specimen transport bag, zipper closure	1000
4929	12 x 15" specimen transport bag, zipper closure	500

Custom printed bags and special sizes are available. Please contact our sales department for information.

SPECIMEN TRANSPORT BAGS - GLUE CLOSURE

The special glue closure of our specimen transport bags offer a liquid-tite tamper evident seal. Simply pull off the liner to expose the glue adhesive and press to close. Item #4923 features a 3 x5" liquid absorbent pad that quickly absorbs aqueous solutions or bodily fluids providing for safe transportation of specimens.

- 1.5 mil gauge polyethylene (PE)
- Attached document compartment
- Glue closure tamper evident seal
- Imprinted with biohazard symbol
- Perfect for the transport of drug testing products



Item #	Description	Unit
4921	6" x 9" specimen transport bag, glue closure	1000
4923	6" x 9" specimen transport bag, glue closure, with absorbent pad	500

BITRAN SERIES S STORAGE BAGS

Bitran® Series S bags are designed for long-term handling and storage of specimens. The double track zipper is easy to close and provides an airtight seal. Produced from 3 mil Saranex®, these bags provide exceptional vapor and odor control.

Item #	Description	Unit
4740	2" x 4", 4 x 250/case	1000
4741	3" x 6", 4 x 250/case	1000
4742	6" x 6", 4 x 250/case	1000
4743	7" x 8", 4 x 250/case	1000
4744	9" x 12", 4 x 250/case	1000
4745	12" x 12", 4 x 250/case	1000
4746	16" x 16", 100/case	100
4747	24" x 24", 100/case	100

