

MICROCENTRIFUGE TUBES

ClickSeal™ Tubes

Research Grade MicroCentrifuge Tubes



Increased volumes of **0.65ml**, **1.7ml** and **2.0ml**, yet fit all popular microcentrifuge rotors. Writing surface on flat cap and side wall for quick sample ID. Unique side-gated mold design eliminates risk of glove or finger puncture caused by sharp bottom-gate protrusions on other tubes. Superior clarity and external graduations ensure accurate, reproducible lab results. Target syringe port is easily pierced for spill-free handling of hazardous substances. With-stands temperatures from +121C to -80C and an RCF of 17,000xg. Positive ClickSeal™ design allows repeated cap closures and eliminates pressure induced popping of caps during incubation. All tubes available in natural polypropylene, in RainbowPacks of five transparent colors, or as packs of a single color. Five bright fluorescent color tubes are offered in **0.65ml**, **1.7ml** and **2.0ml** sizes. Tubes are also bagged and radiation sterilized in small RadPacks™. **Available in the following colors: Natural, Red, Yellow, Green, Violet.**



Rainbow Packs (5 assorted colors)



Sterile Rad Packs

Graduated 0.65ml ClickSeal™ MicroCentrifuge Tubes

Catalog #	Description	Packaging	Qty/Case
CN0650-GT	Graduated 0.65ml ClickSeal MicroCentrifuge Tubes w/Caps Natural color	10 Packs of 1,000	10,000 tubes per case
CN0650-GTS	Sterilized, Graduated 0.65ml ClickSeal MicroCentrifuge Tubes w/Caps, Natural color	10 Bags of 50 Tubes/Pack, 10 Packs of 500 Tubes	5,000 tubes per case
CA0650-GT	Graduated 0.65ml ClickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	10 Packs of 1,000 Tubes	10,000 tubes per case
CA0650-GTS	Sterilized, Graduated 0.65ml ClickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	10 Bags of 50 Tubes/Pack, 10 Packs of 500 Tubes	5,000 tubes per case

CR0650-GT	Red, Graduated 0.65ml ClickSeal Tubes	10 Packs of 1,000 Tubes	10,000 tubes per case
CY0650-GT	Yellow, Graduated 0.65ml ClickSeal Tubes	10 Packs of 1,000 Tubes	10,000 tubes per case
CG0650-GT	Green, Graduated 0.65ml ClickSeal Tubes	10 Packs of 1,000 Tubes	10,000 tubes per case
CB0650-GT	Blue, Graduated 0.65ml ClickSeal Tubes	10 Packs of 1,000 Tubes	10,000 tubes per case
CV0650-GT	Violet, Graduated 0.65ml ClickSeal Tubes	10 Packs of 1,000 Tubes	10,000 tubes per case
CN0640-GT	Graduated 0.65ml ClickSeal MicroCentrifuge Tubes w/Caps Natural color	5 Packs of 1,000	5,000 tubes per case
CN0640-GTS	Sterilized, Graduated 0.65ml ClickSeal MicroCentrifuge Tubes w/Caps, Natural color	10 Bags of 50 Tubes/Pack, 5 Packs of 500 Tubes	2,500 tubes per case
CA0640-GT	Graduated 0.65ml ClickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	5 Packs of 1,000 Tubes	5,000 tubes per case
CA0640-GTS	Sterilized, Graduated 0.65ml ClickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	10 Bags of 50 Tubes/Pack, 5 Packs of 500 Tubes	2,500 tubes per case
CR0640-GT	Red, Graduated 0.65ml ClickSeal Tubes	5 Packs of 1,000 Tubes	5,000 tubes per case
CY0640-GT	Yellow, Graduated 0.65ml ClickSeal Tubes	5 Packs of 1,000 Tubes	5,000 tubes per case
CG0640-GT	Green, Graduated 0.65ml ClickSeal Tubes	5 Packs of 1,000 Tubes	5,000 tubes per case
CB0640-GT	Blue, Graduated 0.65ml ClickSeal Tubes	5 Packs of 1,000 Tubes	5,000 tubes per case
CV0640-GT	Violet, Graduated 0.65ml ClickSeal Tubes	5 Packs of 1,000 Tubes	5,000 tubes per case

Graduated 1.7ml ClickSeal™ MicroCentrifuge Tubes

Catalog #	Product Description	Packaging	Qty/Case
CN1700-GT	Graduated 1.7ml ClickSeal MicroCentrifuge Tubes w/Caps, Natural color	10 Packs of 500 Tubes	5,000 tubes per case
CN1700-GTS	Sterilized, Graduated 1.7ml ClickSeal MicroCentrifuge Tubes w/Caps, Natural color	5 Bags of 50 Tubes/Pack, 10 Packs of 250 Tubes	2,500 tubes per case
CA1700-GT	Graduated 1.7ml ClickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	10 Packs of 500 Tubes	5,000 tubes per case
CA1700-GTS	Sterilized, Graduated 1.7ml ClickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	5 Bags of 50 Tubes/Pack, 10 Packs of 250 Tubes	2,500 tubes per case
CR1700-GT	Red, Graduated 1.7ml ClickSeal Tubes	10 Packs of 500 Tubes	5,000 tubes per case
CY1700-GT	Yellow, Graduated 1.7ml ClickSeal Tubes	10 Packs of 500 Tubes	5,000 tubes per case
CG1700-GT	Green, Graduated 1.7ml ClickSeal Tubes	10 Packs of 500 Tubes	5,000 tubes per case
CB1700-GT	Blue, Graduated 1.7ml ClickSeal Tubes	10 Packs of 500 Tubes	5,000 tubes per case
CV1700-GT	Violet, Graduated 1.7ml ClickSeal Tubes	10 Packs of 500 Tubes	5,000 tubes per case
CN1710-GT	Natural, Graduated 1.7ml ClickSeal Tubes	1 Bag of 1000Tubes/Pack, 10 Packs of 1000 Tubes	10,000 tubes per case
CA1710-GT	Graduated 1.7ml ClickSeal Tubes in RainbowPack	1 Bag of 1000Tubes/Pack, 10 Packs of 1000 Tubes	10,000 tubes per case

Graduated 2.0ml ClickSeal™ MicroCentrifuge Tubes

Catalog #	Product Description	Packaging	Qty/Case
CN2000-GT	Graduated 2.0ml ClickSeal MicroCentrifuge Tubes w/Caps, Natural color	10 Packs of 400 Tubes	4,000 tubes per case
CN2000-GTS	Sterilized, Graduated 2.0ml ClickSeal MicroCentrifuge Tubes w/Caps, Natural color	4 Bags of 50/Pack, 10 Packs of 200 Tubes	2,000 per case
CA2000-GT	Graduated 2.0ml ClickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	10 Packs of 400 Tubes	4,000 tubes per case
CA2000-GTS	Sterilized, Graduated 2.0ml ClickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	4 Bags of 50/Pack, 10 Packs of 200 Tubes	2,000 tubes per case
CR2000-GT	Red, Graduated 2.0ml ClickSeal Tubes	10 Packs of 400 Tubes	4,000 tubes per case
CY2000-GT	Yellow, Graduated 2.0ml ClickSeal Tubes	10 Packs of 400 Tubes	4,000 tubes per case
CG2000-GT	Green, Graduated 2.0ml ClickSeal Tubes	10 Packs of 400 Tubes	4,000 tubes per case
CB2000-GT	Blue, Graduated 2.0ml ClickSeal Tubes	10 Packs of 400 Tubes	4,000 tubes per case
CV2000-GT	Violet, Graduated 2.0ml ClickSeal Tubes	10 Packs of 400 Tubes	4,000 tubes per case

SlickSeal™ Tubes



Pre-Lubricated MicroCentrifuge Tubes

Pre-lubricated microcentrifuge tubes simplify research procedures requiring a low adhesion surface. Ensures optimal sample yield during micro-volume sequencing procedures. Ideal for custom gene libraries, protein work, cDNA synthesis, and nucleic acid amplification. Pre-lubricated SlickSeal tubes incorporate standard ClickSeal design features and packaging options, tubes are offered in **0.65ml, 1.7ml, and 2.0ml** sizes. **Available in the following colors: Natural, Red, Yellow, Green, Violet.**

Pre-Lubricated Graduated 0.65ml SlickSeal™ MicroCentrifuge Tubes

Catalog #	Description	Packaging	Qty/Case
CN065S-GT	Graduated 0.65ml SlickSeal MicroCentrifuge Tubes w/Caps Natural color	10 Packs of 500 Tubes	5,000 tubes per case
CN065S-GTS	Sterilized, Graduated 0.65ml SlickSeal MicroCentrifuge Tubes w/Caps, Natural color	10 Bags of 50 Tubes/Pack, 10 Packs of 250 Tubes	2,500 tubes per case

Pre-Lubricated Graduated 1.7ml SlickSeal™ MicroCentrifuge Tubes

Catalog #	Product Description	Packaging	Retail
CN170S-GT	Graduated 1.7ml SlickSeal MicroCentrifuge Tubes w/Caps, Natural color	10 Packs of 250 Tubes	2,500 tubes per case
CN170S-GTS	Sterilized, Graduated 1.7ml SlickSeal MicroCentrifuge Tubes w/Caps, Natural color	5 Bags of 50 Tubes/Pack, 5 Packs of 250 Tubes	1,250 tubes per case
CA170S-GT	Graduated 1.7ml SlickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	10 Packs of 250 Tubes	2,500 tubes per case
CA170S-GTS	Sterilized, Graduated 1.7ml SlickSeal MicroCentrifuge Tubes w/Caps, RainbowPack	5 Bags of 50 Tubes/Pack, 5 Packs of 250 Tubes	1,250 tubes per case
CR170S-GT	Red, Graduated 1.7ml SlickSeal Tubes	10 Packs of 250 Tubes	2,500 tubes per case
CY170S-GT	Yellow, Graduated 1.7ml SlickSeal Tubes	10 Packs of 250 Tubes	2,500 tubes per case
CG170S-GT	Green, Graduated 1.7ml SlickSeal Tubes	10 Packs of 250 Tubes	2,500 tubes per case
CB170S-GT	Blue, Graduated 1.7ml SlickSeal Tubes	10 Packs of 250 Tubes	2,500 tubes per case
CV170S-GT	Violet, Graduated 1.7ml SlickSeal Tubes	10 Packs of 250 Tubes	2,500 tubes per case
CN171S-GT	Graduated 1.7ml SlickSeal MicroCentrifuge Tubes w/Caps, Natural color, bulk	1 Bags of 500 Tubes/Pack, 10 Bags of 500 Tubes	5,000 tubes per case
CA171S-GT	Graduated 1.7ml SlickSeal MicroCentrifuge Tubes w/Caps, RainbowPack, bulk	1 Bags of 500 Tubes/Pack, 10 Bags of 500 Tubes	5,000 tubes per case

Pre-Lubricated Graduated 2.0ml SlickSeal™ MicroCentrifuge Tubes

Catalog #	Product Description	Packaging	Retail
CN200S-GT	Graduated 2.0ml SlickSeal MicroCentrifuge Tubes w/Caps, Natural Color.	10 Packs of 200 Tubes	2,000 tubes per case
CN200S-GTS	Sterilized, Graduated 2.0ml SlickSeal MicroCentrifuge Tubes w/Caps, Natural Color.	4 Bags of 50/Pack, 20 Packs of 50 Tubes	1,000 tubes per case

CAESA-LAB INC.

Tel: (450) 441-5408 / Fax: (450) 441-9167

www.caesalabinc.com

e-mail: info@caesalabinc.com

