

M920

StainTray™ SLIDE STAINING SYSTEM

Made of ABS Plastic

Another user friendly approach to immuno-histochemistry staining. This tray is also suitable not only for routine staining requiring a humid chamber but is also ideal for Hematology, Cytology and Microbiology laboratories. Manipulation is made safe and easy by using only one hand.

The StainTray™ has a black base made of tough ABS plastic withstanding a wide range of chemicals*. It will accept up to 20 slides on four plastic rails covered with a polymer strip to perfectly hold slides even if tray is held at an angle. When humidity is needed, wells between rails will hold water securely without splashing. Rails are raised not only to avoid water touching the slides but to make them more easily retrieveable. The base will also hold excess stain solution dripping from the slides. Four rubber feet ensure greater base stability. Units are stackable for space saving purposes.

Two covers are available:

- A clear one allowing for visual examination. Made of PETG with a temperature range of -20 °C to 60 °C.
- A black lid for fluorescent work. Made of ABS with a temperature range of -80 °C to 80 °C.

Dimensions with cover: 38 x 24 x 4.5 cm H. (15 x 9 ³/₈ x 1 ³/₄ in. H.)



Unit also available with black cover.

* Avoid chlorinated hydrocarbons.

| Cat.# | Description | Capacity | Qty/Cs |
|--------|-------------------------|-----------|--------|
| M918-1 | Base with Clear Lid | 10 slides | 1 |
| M918-2 | Base with Black Lid | 10 slides | 1 |
| M919-1 | Clear Lid only for M918 | — | 1 |
| M919-2 | Black Lid only for M918 | — | 1 |
| M920-1 | Base with Clear Lid | 20 slides | 1 |
| M920-2 | Base with Black Lid | 20 slides | 1 |
| M921-1 | Clear Lid only for M920 | — | 1 |
| M921-2 | Black Lid only for M920 | — | 1 |



AVAILABLE IN NOVEMBER 2006

A 10-slide StainTray™ offering the same great features and benefits



Plastic rails are covered with a polymer strip to perfectly hold slides even if tray is held at an angle.



Wells between rails can be used to hold water for techniques needing a humid environment.



Drain plug can be removed for easier emptying of SlideTray.