

STORAGE BOXES

T514 StoreBox™ STORAGE BOXES

Made of high impact polystyrene

Color your world with a wide variety of economical storage boxes for tubes from 1.2 ml to 5 ml.


These storage boxes are designed to be used at temperatures between -90 °C and +80 °C. Different models are available to accommodate either 25, 81 or 100 sample tubes.

A transparent cover allows the user to see the contents of the box, and is keyed to the base to prevent misalignment. Printed with a series of squares (numbered from 1 to 25, 1 to 81, or 1 to 100), the surface accepts writing with markers, facilitating inventory control.

A unique color coding system uses colored plastic grids to separate the cover from the base on the 25 and 81-place boxes. Those made to accept 100 tubes (serie 2100) have a colored base instead of a grid. Removal of vials facilitated by an innovative vial picker supplied with each storage box. A choice of four popular pastel colors is available. Not autoclavable.




Series 225, 281 & 581 Storage Boxes



FEATURES AND BENEFITS

- 1 Cover has numbered squares for easy sample identification
- 2 Two corners of cover and base are cut to prevent misalignment
- 3 Writing surface for identifying base and/or cover
- 4 Vials readily visible through transparent cover
- 5 Air vents minimizing condensation
- 6 Drain holes under base
- 7 Stackable
- 8 Four pastel colors available for better color-coding



ALL STOREBOX™ STORAGE BOXES ARE EASILY STACKABLE

StoreBox™ Selection Guide

Please note that tubes and caps are sold separately.



	-1T	-2T	-4T	-4AT	-5T	-1ATPR	-2TPR	-2ATPR	-3ATPR	-4TPR	-4ATPR	-5TPR	-5ATPR
Tube Capacity (ml)	1.2	2	4	4	5	1.2	2	2	3	4	4	5	5
25-place Storage Box Series T514-225	•	•				•	•	•					
81-place Storage Box Series T514-281			•	•		•	•	•					
81-place Storage Box Series T514-581				*	*				•	•	•	•	•
100-place Storage Box Series T514-2100	•	•											

Series 225: 76 mm x 76 mm x 52 mm H. (3 x 3 x 2 1/16 in. H.)				
Cat. No.	For sample tubes	Color of grid	Qty/Pk	Qty/Cs
T514-225B	1 to 2 ml	Blue	8	48
T514-225G	1 to 2 ml	Green	8	48
T514-225P	1 to 2 ml	Pink	8	48
T514-225Y	1 to 2 ml	Yellow	8	48
Series 281: 133 mm x 133 mm x 52 mm H. (5 1/4 x 5 1/4 x 2 1/16 in. H.)				
Cat. No.	For sample tubes	Color of grid	Qty/Pk	Qty/Cs
T514-281B	1 to 2 ml	Blue	4	24
T514-281G	1 to 2 ml	Green	4	24
T514-281P	1 to 2 ml	Pink	4	24
T514-281Y	1 to 2 ml	Yellow	4	24
Series 581: 133 mm x 133 mm x 95 mm H. (5 1/4 x 5 1/4 x 3 3/4 in. H.)				
Cat. No.	For sample tubes	Color of grid	Qty/Pk	Qty/Cs
T514-581B	3 to 5 ml	Blue	5	10
T514-581G	3 to 5 ml	Green	5	10
T514-581P	3 to 5 ml	Pink	5	10
T514-581Y	3 to 5 ml	Yellow	5	10
Series 2100: 133 mm x 133 mm x 52 mm H. (5 1/4 x 5 1/4 x 2 1/16 in. H.)				
Cat. No.	For sample tubes	Color of base	Qty/Pk	Qty/Cs
T514-2100B	1 to 2 ml	Blue	4	24
T514-2100G	1 to 2 ml	Green	4	24
T514-2100P	1 to 2 ml	Pink	4	24
T514-2100Y	1 to 2 ml	Yellow	4	24
Cat. No.	Description		Qty/Cs	
T314S	Metal stand		1	

* For these tubes, see T314-481 polycarbonate boxes on pages 64-65.

Removal of vials facilitated by an innovative vial picker supplied with each StoreBox™.



T314S Metal Stand allows a storage box to be kept at an angle for easier access to tubes.