


Histoplex™

Tissue cassettes with micromesh chamber



Effective for the embedding, processing and storage of both standard tissue specimens and tissue fragments. The innovative micromesh chamber securely retains fragments leaving ample room available for additional larger specimens. Full range of colors available for identification.

s a f e t y , q u a l i t y

& r e l i a b i l i t y

HISTOPLEX

Features:

Durable acetyl polymer resists the chemical action of histological agents

Flow through slots maximize fluid exchange


Fine micromesh chamber, .23mm, ideal for tissue fragments, and effective for fluid exchange, no need for the additional use of filters

Textured, 45° angled writing surface, suitable for labeling devices

Clam-shell design can be snapped closed repeatedly and reopens with ease

Non-cytotoxic, range of 12 non-metallic colors

Conveniently packaged, 500 cassettes per dispenser box, 3 boxes per case



TC-WHITE

Product Codes for Available Colors:

TC-AQUA	TC-BLUE	TC-YELLOW	TC-PINK	TC-ORANGE	TC-GREEN
					
TC-TAN	TC-RED	TC-LILAC	TC-GREY	TC-DRKGREEN	TC-PEACH
					

*White also available

CAESA-LAB INC

TELEPHONE : (450) 441-5408

FAX : (450) 441-9167

E-MAIL : info@caesalabinc.com

www.caesalabinc.com